

temperature of the probe, and 3°C below the melting temperature of the probe] of 37°C
for 16 hours in 5X SSC, 5X Denhardt solution, 25% formamide, and 100 µg/ ml
denatured salmon sperm DNA, wherein said probe comprises an HIV-2 nucleic acid
molecule, which hybridizes to HIV-2ROD genomic DNA under hybridization conditions
[selected from the group consisting of 42°C below the melting temperature of the
probe, 20°C below the melting temperature of the probe, and 3°C below the melting
temperature of the probe] of 37°C for 16 hours in 5X SSC, 5X Denhardt solution, 25%
formamide, and 100 µg/ ml denatured salmon sperm DNA, with washes in 2X SSC,
0.1% SDS at 25°C; 1X SSC, 0.1% SDS at 60°C; or 0.1X SSC, 0.1% SDS at 60°C;

b) washing the resulting hybrid under conditions selected from the group
consisting of [42°C below the melting temperature of the probe, 20°C below the melting
temperature of the probe, and 3°C below the melting temperature of the probe] 2X
SSC, 0.1% SDS at 25°C; 1X SSC, 0.1% SDS at 60°C; and 0.1X SSC, 0.1% SDS at
60°C; and

c) detecting said hybrid.

92. (Amended) A method of detecting HIV-2 retrovirus nucleic acid in a
biological sample, said method comprising:

a) contacting said sample with an HIV-2 specific probe under
hybridization conditions [selected from the group consisting of hybridization conditions
of 42°C below the melting temperature of the probe, 20°C below the melting
temperature of the probe, and 3°C below the melting temperature of the probe] of 37°C

Ser. No. 08/470,489

RECEIVED
AUG - 5 1993
TC 1600 MAIL ROOM

for 16 hours in 5X SSC, 5X Denhardt solution, 25% formamide, and 100 μ g/ml

denatured salmon sperm DNA,

wherein said probe comprises an HIV-2 nucleic acid molecule obtained from nucleotides 1-380 of the U3/R region of HIV-2, nucleotides 1-1566 of the *gag* gene of HIV-2, nucleotides 1114-1524 of the *gag* gene, nucleotides 1-405 of the *gag* gene, nucleotides 406-1155 of the *gag* gene, or nucleotides 1-2673 of the *env* gene of HIV-2, and

wherein said probe hybridizes to HIV-2ROD genomic DNA under hybridization conditions [selected from the group consisting of 42°C below the melting temperature of the probe, 20°C below the melting temperature of the probe, and 3°C below the melting temperature of the probe] of 37°C for 16 hours in 5X SSC, 5X Denhardt solution, 25% formamide, and 100 μ g/ ml denatured salmon sperm DNA, with washes in 2X SSC, 0.1% SDS at 25°C; 1X SSC, 0.1% SDS at 60°C; or 0.1X SSC, 0.1% SDS at 60°C;

b) washing the resulting hybrid under conditions selected from the group consisting of [42°C below the melting temperature of the probe, 20°C below the melting temperature of the probe, and 3°C below the melting temperature of the probe] 2X SSC, 0.1% SDS at 25°C; 1X SSC, 0.1% SDS at 60°C; and 0.1X SSC, 0.1% SDS at 60°C; and

c) detecting said hybrid.

99. (Amended) A method of producing an HIV-2 specific hybridization probe for HIV-2 retrovirus nucleic acid, said method comprising:

a) preparing a nucleic acid insert, which hybridizes to HIV-2ROD genomic DNA under hybridization conditions [selected from the group consisting of 42°C below the melting temperature of the insert, 20°C below the melting temperature of the insert, and 3°C below the melting temperature of the insert] of 37°C for 16 hours in 5X SSC, 5X Denhardt solution, 25% formamide, and 100 µg/ ml denatured salmon sperm DNA, with washes in 2X SSC, 0.1% SDS at 25°C; 1X SSC, 0.1% SDS at 60°C; or 0.1X SSC, 0.1% SDS at 60°C;

- b) introducing the insert into a recombinant cloning vector;
- c) introducing said vector into a competent cellular host; and
- d) recovering the DNA recombinants.

101. (Amended) A method of producing an HIV-2 specific hybridization probe for HIV-2 retrovirus nucleic acid, said method comprising:

a) preparing a nucleic acid insert, wherein said insert is obtained from nucleotides 1-380 of the U3/R region of HIV-2, nucleotides 1-1566 of the *gag* gene of HIV-2, nucleotides 1114-1524 of the *gag* gene, nucleotides 1-405 of the *gag* gene, nucleotides 406-1155 of the *gag* gene, or nucleotides 1-2673 of the *env* gene of HIV-2, and wherein said insert hybridizes to HIV-2ROD genomic DNA under hybridization conditions [selected from the group consisting of 42°C below the melting temperature of the insert, 20°C below the melting temperature of the insert, and 3°C below the melting